cTn I Reag

058615/01-1506

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1. Tobacman LS. Thin Filament-Mediated Regulation of Cardiac Contraction. *Annu Rev Physiol*. 58:447-481; 1996.
2. Wilkinson JM, Grand JA. Comparison of Amino Acid Sequence of TnI from Different Striated Muscles. *Nature*. 271:31-35; 1978.

3. Bodor GS. Cardiac TnI: A Highly Specific Biochemical Marker for Myocardial Infarction. *J Clin Immunoassay*. 17(1):40-44; 1994.
4. Vallins WJ et al. Molecular Cloning of Human cTnI Using Polymerase Chain Reaction. *FEBS Letts*. 270(1,2):57-61; 1990.
5. Mair J et al. Equivalent Early Sensitivities of Myoglobin, CK-MB Mass, CK Isoform Ratios and Cardiac Troponins I and T for AMI. *Clin Chem*. 41(9):1266-1272; 1995.
6. Lerue C et al. Cardiac-Specific Immunoassayometric Assay of TnI in the Early Phase of AMI. *Clin Chem*. 39(6):972-979; 1993.
7. Adams JE et al. Cardiac Troponin I. A Marker with High Specificity for Cardiac Injury. *Circulation*. 88(1):101-106; 1993.
8. Apple FS et al. Cardiac Troponin, CK-MB and Myoglobin for the Early Detection of AMI and Monitoring of Reperfusion Following Thrombolytic Therapy. *Clin Chim Acta*. 237:59-66; 1995.
9. Jaffe AS. Elevations in cardiac troponin measurements: false-positive. *Cardio Tox* 2001; 1:87-92.
10. Galvani M, Ottani F, Ferrini D, et al. Prognostic Influence of Elevated Values of Cardiac Troponin I in Patients with Unstable Angina. *Circulation*. 95:2053-2059; 1997.
11. Antman EM, Tanasijevic MJ, Thompson B, et al. Cardiac-Specific Troponin I Levels to Predict the Risk of Mortality in Patients with Acute Coronary Syndromes. *N Eng J Med*. 335:1342-1349; 1996.
12. Keffer JH. The Cardiac Profile and Proposed Practice Guideline for Acute Ischemic Heart Disease. *Am J Clin Pathol*. 107:398-409; 1997.
13. Missov E, Calzolari C, Davy JM, Leclercq F, et al. Cardiac troponin I in patients with hematologic malignancies. *Coron Artery Dis*. 8:537-541; 1997.
14. *Markers in Cardiology: Current and Future Clinical Applications*. Edited by Adams JE, Apple FS, Jaffe AS, Wu AHB (ISBN0-87993-472-7). Futura Publishing, Armonk NY 10504; 2001.
15. Lindahl B, Diderholm E, Langerquist B, et al, for the FRISC II Investigators. Mechanisms behind the prognostic value of troponin T in unstable coronary artery disease: FRISC II substudy. *J Am Coll Cardiol*. 2001;38: 1497-1498.
16. Venge P, Langerquist B, Diderholm E, et al, on behalf of the FRISC II Study Group. Clinical Performance of three cardiac troponin assays in patients with unstable coronary artery disease (FRISC II substudy). *Am J Cardiol*. 2002;89: 1035-1041.
17. Ottani F, Galvani M, Nicolini A, et al. Elevated cardiac troponin levels predict the risk of adverse outcome in patients with acute coronary syndrome. *Am Heart J*. 2000;40:917-927.
18. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin in patients with non-ST elevation acute coronary syndromes: a meta analysis. *J Am Coll Cardiol*. 2001;38:478-485.
19. Summers M et al. Luminogenic Reagent Using 3-chloro 4-hydroxy acetanilide to Enhance Peroxidase/Luminol Chemiluminescence. *Clin Chem*. 41:S73; 1995.
20. NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition*. NCCLS document H3-A5 (ISBN 1-56238-515-1). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, 2003.
21. NCCLS. *Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Fifth Edition*. NCCLS document H4-A5 [ISBN 1-56238-538-0]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
22. CLSI. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guidelines-Second Edition*. CLSI Document C28-A2 (ISBN 1-56238-406-6); CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087, 2000.
23. CLSI. *Assessment of the Clinical Accuracy of Laboratory Tests Using Receiver Operator Characteristics (ROC) Plots; Approved Edition*; CLSI Document GP10-A (ISBN 1-56238-285-3); CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087, 2001.
24. Joint European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction Redefined3/4A Consensus Document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. *J Am Coll Cardiol*. 36:959-969; 2000.
25. Braunwald E, Antman EM, Beasley JW, Califf RM, Chaitlin MD, Hochman JS, Jones RH, Kereiakes D, Kupersmith J, Levin TN, Pepine CJ, Schaeffer JW, Smith EE, Steward DE, Teroux P. ACC/AHA guidelines for the management of patients with unstable angina and non-ST-segment elevation myocardial infarction: Executive summary and recommendations, A report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on the

- Management of Patients with Unstable Angina). *Circulation* 2000; 102(10):1193-1209.
26. Braunwald E, Antman EM, Beasley JW, Califf RM, Cheitlin MD, Hochman JS, Jones RH, Kereiakes D, Kupersmith J, Levin TN, Pepine CJ, Schaeffer JW, Smith EE, Steward DE, Teroux P. ACC/AHA guidelines for the management of patients with unstable angina and non-ST-segment elevation myocardial infarction: a report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines (Committee on the Management of Patients with Unstable Angina), 2002.
 27. Wu AH, Apple FS, Gibler WB, Jesse RL, Warshaw MM, Valdes R Jr. National Academy of Clinical Biochemistry Standards of Laboratory Practice: Recommendations for the Use of Cardiac Markers in Coronary Artery Diseases. *Clinical Chemistry*.45:1104-1121; 1999.
 28. Panteghini M, Apple FS, Christenson RH, Dati F, Mair J, Wu AH. Use of Biochemical Markers in Acute Coronary Syndromes. IFCC Scientific Division, Committee on Standardisation of Markers of Cardiac Damage. *Clin Chem Lab Med.* 37:687-693; 1999.
 29. NACB Laboratory Medicine Practice Guidelines: Characteristics & Utilization of Biochemical Markers in ACS and Heart Failure CLINICAL: ACUTE CORONARY SYNDROMES Chapter 1 Version Date April 8, 2004.
 30. World Health Organization. Report of the Joint International Society and Federation of Cardiology/WHO Task Force on Standardization of Clinical Nomenclature. Nomenclature and Criteria for Diagnosis of Ischemic Heart Disease. *Circulation.* 59:607-609; 1979.
 31. NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline.* NCCLS document EP17-A (ISBN 1-56238-551-8). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087, 2004.
 32. NCCLS. *Evaluation of Precision performance of Quantitative Measurement Methods; Approved Guideline-Second Edition.* NCCLS document EP5-A2 (ISBN 1-56238-542-9). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA 2004
 33. NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.* NCCLS document EP9-A2 (ISBN 1-56238-472-4). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2002.
 34. Deming WE. *Statistical Adjustment of Data.* New York, NY: John Wiley and Sons; 1943.

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